

# Product Information

## Cyanine 555-dUTP, Lyophilized powder

**Catalog Number** 40064

**Unit Size** 25 nmole

**Molecular Weight** ~1131

**Color and Form** Red Solid

**Absorption/Emission** 550/570 nm

### Storage and Handling

Store Cyanine 555-dUTP at -20°C or -70°C, protected from light. When stored as recommended, product is stable for at least 6 months from date of receipt.

### Product Application

Fluorescent conjugates of dUTP can be used for detection of apoptotic cells by direct fluorescence TUNEL labeling of DNA strand breaks in cells<sup>1</sup>, and can be used in place of dTTP in standard DNA labeling and synthesis protocols to generate fluorescent dsDNA and oligonucleotide probes.

### References

1. Gold et al. (1994). *Lab Invest.* 71 (2):219-25.
2. Slupphaug et al. (1993). *Anal Biochem.* 211 (1):164-9.
3. Hogrefe et al. (2002). *PNAS* 99 (2): 596-601.

Please visit our website at [www.biotium.com](http://www.biotium.com) to view our full selection of CF™ dye bioconjugates, including antibodies, antibody labeling kits, phalloidin, annexin V and a-bungarotoxin, as well as classic fluorescent nucleic acid dyes and antifade fluorescence mounting media.