



FlexiFuge™

Put A New Spin on Things With the Argos Brand Flexifuge!

FlexiFuge is a versatile, personal micro-centrifuge for quick, low speed centrifugation. Whisper quiet motor offers two speeds (2000g or 5000g / 6,000 or 9,200 RPM) depending on application. Flexifuge can be operated in continuous mode, pulsed for a short burst of rapid centrifugation, or timed to spin from 1-15 minutes. Regardless of setting, rotors stop once the lid release button is pushed. A ring of LED lights illuminate when the unit is active and can be viewed through the transparent lid. Color can be changed (red, green, blue) or turned off by a switch under the unit.

Flexifuge is supplied with an 8-Place fixed angle rotor for 1.5mL tubes, 4-Place rotor for 5mL tubes, rotor for (4), 8 x.2mL PCR tube strips and 8 adapters for .5mL tubes.

LED lights glow Red, Green or Blue

or can be turned off by switch on bottom of unit



FlexiFuge by Argos

The features you've come to expect

- Whisper quiet motor
- Two speeds - 6,000 or 9,200 RPM's
- Operates in continuous mode, pulse or timer (1-15 min.)
- Accommodates **NEW 5mL rotor!**



FlexiFuge can be operated in continuous mode or pulsed for rapid centrifugation



Includes 8-Place fixed angle rotor for 1.5mL tubes, 4-Place rotor for 5mL tubes, rotor for (4) 8 x 0.2mL PCR tube strips and 8 adapters for 0.5mL tubes



Timer can be set for 1-15 minutes



Argos' Exclusive 5mL Tube!

A Perfect Compliment to the FlexiFuge...

The 5mL polypropylene tube is the ideal alternative to 15mL conical tubes when using volumes less than 5mL. Its smaller size means you no longer need to insert your pipettor barrel inside the tube. Withstands boiling, autoclaving (120°C for 15 minutes at 15 psi), and is suitable for freezing to -90°C. Sterile tubes available upon request.

Cat. No.	Description		Qty.
C1000	Flexifuge Centrifuge 120V	C1000	1/EA
C0300	1.5mL Rotor	C0300	1/EA
C0301	PCR Strip Rotor	C0301	1/EA
C0302	5mL Rotor	C0302	1/EA
C0303	Adapters for .5mL Tubes, 8 each	C0303	1/PK
T2076	5mL Conical Prep Tube, Natural	T2076	250/PK
T2077	5mL Conical Prep Tube, Blue	T2077	250/PK
T2078	5mL Conical Prep Tube, Yellow	T2078	250/PK
T2079	5mL Conical Prep Tube, Green	T2079	250/PK
T2080	5mL Conical Prep Tube, Rose	T2080	250/PK
T2081	5mL Conical Prep Tube, Black	T2081	250/PK



Ph: 800-886-8675 • Fx: 847-622-0454
www.argos-tech.com

Distributed by: |

REV0809 LIT.NO. 7000

