



**New Grid Box Design!**



**PolarSafe 81 & 100 Place Gridded Cardboard Box**

Argos now offers our most popular box sizes with a white, water resistant, unique paper lining on the inside box top and bottom; as well as the cell dividers. This water barrier prevents water from penetrating the box, which can cause decomposition and mold.

Each box also contains a numerical grid that allows the researcher to replace the tube exactly where it was removed. A great tool to help keep your samples organized.

**PolarSafe PCR Cardboard Box**

The newest addition to the Argos cryostorage line features a 5 1/4" x 5 1/4" x 1" cardboard PCR box. The box will accommodate a 196-place cardboard divider; included. The box is designed to hold single PCR tubes or PCR strips. The 1" height prevents the tubes from dropping down into the box as often occurs when using a standard 2" box.

**Cryo/Freezer Boxes with Cell Dividers**

Dimensions L x W x H	Number of Cells	Catalog Number
3 x 3 x 2"	25	3025A
5 1/4 x 5 1/4 x 2"	81 (Box w/grid)	R3026A
5 1/4 x 5 1/4 x 2"	100 (Box w/grid)	R3027A
5 1/4 x 5 1/4 x 1"	196	R3012A

Sold individually

**Cell Dividers only for 15mL Tubes**

Dimensions L x W x H	No. of Spaces	Cell Opening	Catalog Number
5 3/8 x 5 3/8 x 2 3/4"	36	22.0mm	R3022

Sold 12 per pack

**Cryo/Freezer Boxes without Cell Dividers**

Dimensions L x W x H	Catalog Number
3 x 3 x 2"	R3013A
5 3/4 x 5 3/4 x 4 7/8"	R3014A
5 1/4 x 5 1/4 x 2"	R3015A
5 1/4 x 5 1/4 x 3"	R3016A

Sold individually

**Cell Dividers only for Standard Tubes**

Dimensions L x W x H	No. of Spaces	Cell Opening	Catalog Number
2 5/8 x 2 5/8 x 1 1/4"	25	11.5mm	R3028
4 7/8 x 4 7/8 x 1 5/16"	25	24.5mm	R3023
4 7/8 x 4 7/8 x 1 5/16"	49	16.7mm	R3020
4 7/8 x 4 7/8 x 1 5/16"	64	14.6mm	R3019
4 7/8 x 4 7/8 x 1 5/16"	81	12.7mm	R3018
4 7/8 x 4 7/8 x 1 5/16"	100	11.5mm	R3017
4 7/8 x 4 7/8 x 1 5/16"	144	9.5mm	R3024
4 7/8 x 4 7/8 x 3/4"	196	8mm	R3011

Sold 12 per pack

**Cell Dividers only for 50mL Tubes**

Dimensions L x W x H	No. of Spaces	Cell Opening	Catalog Number
5 3/8 x 5 3/8 x 2 3/4"	16	33.5mm	R3021

Sold 12 per pack