

BRAND® Transferpette® electronic

Single & Multichannel Pipettes

German manufactured; unprecedented quality and value in an electronic pipette



World's first independently certified "Ergonomics Approved" single and multichannel pipettes!



Intuitive Interface

User-friendly design allows easy operation, even for new lab techs.

Adjustable Finger Rest


The adjustable finger rest adapts the instrument to your most comfortable hand position.

Rechargeable Battery

Over 4000 pipetting cycles between charges means well over a full day of uninterrupted pipetting.

Multichannel V-rings

V-rings on tip cones of multichannel models improve tip fit and reduce mounting and ejection forces.



Easy Maintenance

Multichannel manifolds disassemble easily for cleaning or service. See page 6 for more details on the manifolds of BRAND'S mechanical and electronic pipettes.

Transferpette® electronic single channel pipettes

Capacity, μ L	without charger		with US plug charger	
	Cat. No.	2017 List Price	Cat. No.	2017 List Price
0.5-10	705339	\$479.00	705319	\$510.00
2-20	705340	479.00	705320	510.00
20-200	705343	479.00	705323	510.00
100-1,000	705346	479.00	705326	510.00
500-5,000 (0.5-5mL)	705347	479.00	705327	510.00

Transferpette®-8 and -12 electronic multichannel pipettes

Capacity, μ L	8-channel		12-channel	
	Cat. No.	2017 List Price	Cat. No.	2017 List Price
0.5-10	705419	\$917.00	705469	\$1,097.00
1-20	705420	917.00	705470	1,097.00
5-100	705423	917.00	706573	1,097.00
10-200	705424	917.00	705474	1,097.00
15-300	705426	917.00	705476	1,097.00

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726

www.brandtech.com