PDS No. 6391xx/ 639160-TRI	PRODUCT DATA SHEET			Page 1 of 1	
Revision 11	Petri Dish, 145 x 20 mm			6	
Revision 11	Greiner Item-No. 6391xx / 639160-TRI			greiner bio-one	
Valid for Item-No.:	639102	639160 (sterile)	639161 (sterile)	639160-TRI (st	erile)

1.	Description / Specification			
1.1	Description	Petri Dish with vents		
		639102: standard		
		639160: CELLSTAR®, physical surface treatment, sterile		
		639161: sterile		
		639160-TRI: CELLSTAR®, physical surface treatment, sterile, triple packed		
1.2	Dimensions	See Customer Drawings		
		Weight: 639102, -161: 42,0 – 46,0 g (total weight)		
		639160, -160-TRI: 21,1-22,5 g (lid)		
		24,0-25,4 g (dish)		
1.3	Volume	Max. volume: 240 ml		
		Working volume: 25 - 27 ml		
		639160, -160-TRI: growth area: 143 cm ²		
1.4	Material / Resin	<u>Dish</u> : PS (Polystyrene), free of heavy metal		
		<u>Lid</u> : PS (Polystyrene), free of heavy metal		
1.5	Colour	<u>Dish</u> : clear		
		Lid: clear		
1.6	Sterilisation	639102: no		
		639160, -161: SAL 10 ⁻³		
		639160-TRI: SAL 10 ⁻⁶		
1.7	Quality Control	Raw Material-Control: physical testing		
		Product-Control: testing of attributive and variable characteristics in		
		accordance with the valid specification		
1.8	Other Information	- For single use only		
		- Elevations for ventilation of culture		

2.	Features			
2.1	Basic features	39102, -161: -		
		639160, -160-TRI: free of detectable DNase/RNase, human DNA and		
		pyrogens. Contents non-cytotoxic		
2.2	Temperature range	-20°C to +60°C		
2.3	Autoclavability	No		
2.4	Centrifugation, max. RCF	N/A		
2.5	Chemical Resistance	See homepage:		
		https://www.gbo.com/en_INT/know-how-services/download-center.html		
2.6	Shelf life	639102, -161: 5 years after month of production		
		639160, -160-TRI: 4 years after month of production		
2.7	Other Information	-		

3.	Packaging	639160	639102, -161	639160-TRI
3.1	Pieces / Bag	5	15	5 (triple packed)
3.2	Pieces / Box	120	120	120
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality	-	Certificate of Quality

4.	Other Information
4.1	639160 / 639160-TRI: Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This
10	31 July 2015	3 August 2015		document may not be reproduced for any
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
03.03.2015	S. Kaelberer	Dr. L. Breth	Dr. A. Ganser	