

PDS No.75025x	PRODUCT DATA SHEET				Page 1 of 1
Revision 03	Sapphire Pipette Tips, 1250 µl, non-sterile				 greiner bio-one
	Greiner Item-No. 75025x				
Valid for Item-No.:	750250	750254	750255	750257	750258

1.	Description / Specification	
1.1	Description	PP Sapphire Pipette Tips, 1250 µl, graduated, thin-walled top of the tips for reliable fit and optimal seal, non-sterile 750250: Standard Pipette Tips, bagged 750254: Standard Pipette Tips, placed in racks 750255: Low Retention Pipette Tips, placed in racks 750257: Standard Pipette Tips, refill unit 750258: Low Retention Pipette Tips, refill unit The Sapphire Low Retention Pipette Tips feature specially optimised surface properties for maximum precision Matching rack, Cat.-No. 970350 (blue insert)
1.2	Dimensions	Pipette tip: length: 98.4 mm outer diameter: 8.9 mm Rack: approx. 126 x 87 x 110 mm (L x W x H)
1.3	Volume	1250 µl
1.4	Material / Resin	Pipette tip / rack: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: natural Rack / lid: transparent; Rack insert: blue
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	Warehouse conditions: +3°C to +40°C (relative humidity 15-60%)
2.3	Autoclavability	Pipette tip / rack: Yes (121°C / 20 min.)
2.4	Centrifugation, max. RCF	-
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	750250	750254, -255, -257. -258
3.1	Pieces / packaging unit	768 (bag)	96 (rack / refill unit)
3.2	Pieces / box	768	768
3.3	Pieces / case	7680	3840
3.4	Lot-No.	Consists of the work order numbers for the manufacturing team	
3.5	Other Information	-	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 3 July 2015	Date 3 July 2015	Date 3 July 2015	
Date 01.06.2015	Name S. Kaelberer	Name Dr. R. Heller	Name Dr. A. Ganser	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.