I.W. Tremont Co., Inc.

Filter & Technical Specialty Papers

18 Utter Avenue - Hawthorne, New Jersey 07506

Tel: 973-427-3800 Fax: 973-427-3778 www.iwtremont.com

Technical Data	Sheet	Material Designation	Grade F	
Material Properties Image: Binderless Image: Organic Binder Image: Double Laminated Summary Image: Acrylic Binder Image: Laminated Image: Hydrophobic				
This is an ultra high efficiency particulate air filter material especially suited to applications requiring a very high degree of efficiency. The absence of organic binder and fungicide makes this grade eminently suitable for chromatography and analytical filtration.				
Binderless media				
Micron rating	Basis Weight	Caliper Thickness	Mean Pore Size	
0.7	45	0.016	2.8	
μm	lbs/3,000 ft	inches - 4 psi	μm	
	TAPPI Method T410	TAPPI Method T411		
DOP Smoke Penetration	Air Flow Resistance	Tensile Strength MD	Tensile Strength CD	
0.001	51	3.0	2.5	
% at 0.3 µm @	$mm H_2 O @$	lbs / inches	lbs / inches	
10.5 ft/minute ASTM Method D-2986	10.5 ft/minute ASTM Method D-2986	TAPPI Method T494	TAPPI Method T494	
ASTM Method D-2960				
Dry Elongation MD	Dry Elongation CD	Frazier Permeability	Gurley Stiffness	
2.5	2.5	-	-	
%	%	$ft^3/min/ft^2$	mg	
TAPPI Method T494	TAPPI Method T494	0.5 in H_2 O W.G.	TAPPI Method T543	
ASTM Method F778-82				
Water Repellency	Ignition Loss	Comments: Widely used in airborne particulate monitoring		
-	Binderless	applications. Materi efficiency.	al demonstrates a 99.999%	
Inches H ₂ O	% Loss			
2 -				

Actual filtration performance, i.e. efficiency and dust holding capacity, will vary depending upon filter design parameters and the normal variation of the media properties consistent with the specification range. We continuously strive to define our products and hence the specifications are subject to change.