

## BagLight® Multilayer® U

Round bottom blender bag

Ref 132 125 (400 mL)



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

**BagLight® Multilayer® U** is a round bottom bag designed for circulating blenders. With reinforced weldings, it is ultra-resistant. Designed for circulating blenders.

### TECHNICAL SPECIFICATIONS

- Blender bag with round bottom
- Special for circulating blenders
- Reinforced weldings
- **Multilayer®**: multicoated reinforced complex
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- Approved for food contact: Regulation (EC) No 1935/2004
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Average shelf-life: 25 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 400 mL
- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

### BagLight® Multilayer® U 400

Ref 132 125

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 190 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 5 kg
- Pack of 25
- Box of 500