


RollBag®

Sampling bag

Ref 145 040 (1300 mL)



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

RollBag® is a bag for sampling and liquid and solid samples homogenization, up to 10 g.

TECHNICAL SPECIFICATIONS

- Sampling bag
- Ideal for the sampling and blending of any liquid or solid sample
- Wire-reinforced sealing, stainless steel
- To sample and blend in the same bag
- Pre-cut opening of the bag
- Marking zone
- Stands upright alone
- Resistant up to 80°C (176°F)
- Approved for food contact: Regulation (EC) No 1935/2004
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Average shelf-life: 25 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 1300 mL
- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

RollBag® 1300

Ref 145 040

- Max sampling volume: 1300 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 290 mm
- Box dimensions: 34 x 28 x 11 cm, weight: 4.2 kg
- Box of 500