


### Manufacturer of Quality Dishes Since 1968



Rim Feature	Code	Nominal Size (mm)	Style	ISO Mark Target		
Stackable	2901	60 x 15	Mono	No		



Dishes/	Dishes/ Box	Box/ Pallet	Box Dimensions			Вох	Вох	
Sleeve			L	W	Н	Cubic	Weight	
20	500	72	12.25 in	12.50 in	12.50 in	1.11 ft³	8.50 lbs	
20	500		31.12 cm	31.75 cm	31.75 cm	0.03 m <sup>3</sup>	3.86 kgs	

# KORD-Valmark

kord-valmark.com

**Bioplast Manufacturing, L.L.C.** 128 Wharton Road - Bristol, PA - 19007

Customer Service: Toll Free: 877.877.9680 E-mail: sales@bioplastmfg.com

#### ce: 77.9680

#### 2901 - 60 x 15 mm Stackable Petri Dish

- Smaller diameter uses less media "Media Saver"
- Ideal for water treatment and pharmaceutical membrane filter tests
- Clear viewing surface
- Clear top and bottom viewing surfaces
- Disposable (single use) polystyrene petri dish for general microbiological use
- Molded from medical-grade virgin polystyrene
- Optically clear
- Radiation sterilized
- Contents sterile if sleeve is unopened
- Ventilation ribs allow for free air circulation and reduce condensation during incubation
- Consistent flatness allows for even media distribution
- Flared lid skirt and squared corners for easy one-hand operation, especially with gloves
- Certificates of Sterilization Processing available upon request
- Traceable by lot and SKU number
- Certificates of Compliance available upon request

\*media not included






Manufactured from virgin-grade polystyrene for optimum clarity and quality



Ventilation ribs allow for free air circulation and reduce condensation during incubation



Manufactured to the highest quality for consistent flatness and even media distribution





## Manufacturer of Quality Dishes Since 1968



Manufactured to the highest quality for consistent flatness – radiation sterilized Ventilation Ribs allow for full air circulation and reduce condensation during incubation


Flared lid Skirt and squared corners for easy one-hand operation, especially with gloves

Nominal	K-V		Rim Feature	Unique Feature	Dishes/	Dishes/	Cases/ Skid
Size (mm)	Code				Sleeve	Case	
60 x 15	2901	Mono	Stackable	Media Saver	20	500	72
100 x 15	2800	Mono	Fully Stackable	Stable Stacking Design	25	500	30
100 x 15	2900	Mono	Stackable	The Industry Standard	25	500	30
100 x 15	2903	Ві	Stackable	Two Compartments	25	500	30
100 x 15	2912	Tri	Stackable	Three Compartments	25	500	30
100 x 15	2913	Quad	Stackable	Four Compartments	25	500	30
100 x 25	2906	Mono	Stackable	Extra Deep Dish	17	340	30



Automation Petri Dishes (Slippable / No Rim Lid Designs)										
Nominal Size (mm)	K-V Code	Style	Rim Feature	Unique Feature	ISO Mark Target	Dishes/ Sleeve	Dishes/ Case	Cases/ Skid		
60 x 15	2916	Mono	Slippable	Media Saver	No	20	500	72		
95 x 15	95 x 15 2914 Mono Slippable		Media Saver	No	25	500	30			
100 x 10 2904 Mono		Slippable	Space Saver	No	30	750	24			
100 x 15	1920	Mono	No Rim	Automated Lid	No	25	480	30		
100 x 15	2910	Mono	Slippable	The Industry Standard	Yes	25	500	30		
100 x 15	2917	Mono	Slippable	Clear Viewing Surfaces	No	25	500	30		
100 x 15	2920	Mono	No Rim	Automation Lid	No	20	500	30		
100 x 15	2911	Ві	Slippable	Two Compartments	Yes	25	500	30		
100 x 20	2905	Mono	Slippable	Deep Dish, 20mm	Yes	20	400	30		
100 x 25	2909	Mono	Slippable	Extra Deep for Seed Germination Tests	No	17	340	30		
85 x 13	2950	Mono	Fully Stackable/Slippable	New 85 x 13 Ultra Plate	Yes	30	600	30		
150 x 15	x 15 2902 Mono No Rim		No Rim	Concentric Bottom Ring Ensures Even Media Depth	No	10	100	40		





**Contact Us!** 

**Bioplast Manufacturing, L.L.C.** 

128 Wharton Road - Bristol, PA - 19007

**Customer Service:** 

Toll Free: 877.877.9680

E-mail: sales@bioplastmfg.com

