
KORD-Valmark

by Bioplast Manufacturing

Manufacturer of Quality Dishes Since 1968



F	Rim eature	Code	Nominal Size (mm)	Style	ISO Mark Target
SI	lippable	2910	100 x 15	Mono	Yes



Dishes/	Dishes/ Box/		Box Dimensions			Вох	Box	
Sleeve	Box	Pallet	L	W	Н	Cubic	Weight	
25	25 500	500 30	19.25 in	15.50 in	16.00 in	2.76 ft³	20.19 lbs	
25			48.90 cm	39.37 cm	40.64 cm	0.08 m³	9.16 kgs	



kord-valmark.com

Bioplast Manufacturing, L.L.C. 128 Wharton Road - Bristol, PA - 19007

Customer Service: Toll Free: 877.877.9680 E-mail: sales@bioplastmfg.com

2910P - 100 x 15 mm Slippable Petri Dish

- Standard slippable petri dish
- Designed for automated filling and streaking systems
- Side Arrows and ISO Mark Target for automation
- Disposable (single use) polystyrene petri dish for general microbiological use
- Molded from medical-grade virgin polystyrene
- Optically clear
- Radiation sterilized
- Sleeve is printed with the words, "Contents Sterile if Sleeve is Unopened"
- Ventilation ribs allow for free air circulation and reduce condensation during incubation
- Consistent flatness allows for even media distribution
- Flared lid skirt and squared corners for easy one-hand operation, especially with gloves
- Certificates of Sterilization Processing available upon request
- Traceable by lot and SKU number
- Certificates of Compliance available upon request


*media not included



Manufactured to the highest quality for consistent flatness and even media distribution



ISO Mark Target for orientation in automated filling and streaking equipment



Specially designed "slippable" lid is ideal for use in automated equipment



Labware Products by Bioplast Manufacturing

STACKABLA

Manufactured to the highest quality for consistent flatness – radiation sterilized

KORD-Valmark

Ventilation Ribs allow for full air circulation and reduce condensation during incubation Flared lid Skirt and squared corners for easy one-hand operation, especially with gloves

Manufacturer of Quality Dishes Since 1968

Stackable Petri Dishes for General Use							
Nominal Size (mm)	K-V Code	Style	Rim Feature	Unique Feature	Dishes/ Sleeve	Dishes/ Case	Cases/ Skid
60 x 15	2901	Mono	Stackable	Media Saver	20	500	72
100 x 15	2800	Mono	Fully Stackable	Stable Stacking Design	25	500	30
100 x 15	2900	Mono	Stackable	The Industry Standard	25	500	30
100 x 15	2903	Bi	Stackable	Two Compartments	25	500	30
100 x 15	2912	Tri	Stackable	Three Compartments	25	500	30
100 x 15	2913	Quad	Stackable	Four Compartments	25	500	30
100 x 25	2906	Mono	Stackable	Extra Deep Dish	17	340	30



Automation Petri Dishes (Slippable / No Rim Lid Designs)								
Nominal Size (mm)	K-V Code	Style	Rim Feature	Unique Feature	ISO Mark Target	Dishes/ Sleeve	Dishes/ Case	Cases/ Skid
60 x 15	2916	Mono	Slippable	Media Saver	No	20	500	72
95 x 15	2914	Mono	Slippable	Media Saver	No	25	500	30
100 x 10	2904	Mono	Slippable	Space Saver	No	30	750	24
100 x 15	1920	Mono	No Rim	Automated Lid	No	25	480	30
100 x 15	2910	Mono	Slippable	The Industry Standard	Yes	25	500	30
100 x 15	2917	Mono	Slippable	Clear Viewing Surfaces	No	25	500	30
100 x 15	2920	Mono	No Rim	Automation Lid	No	20	500	30
100 x 15	2911	Bi	Slippable	Two Compartments	Yes	25	500	30
100 x 20	2905	Mono	Slippable	Deep Dish, 20mm	Yes	20	400	30
100 x 25	2909	Mono	Slippable	Extra Deep for Seed Germination Tests	No	17	340	30
85 x 13	2950	Mono	Fully Stackable/Slippable	New 85 x 13 Ultra Plate	Yes	30	600	30
150 x 15	2902	Mono	No Rim	Concentric Bottom Ring Ensures Even Media Depth	No	10	100	40



	Contact Us! Bioplast Manufacturing, L.L.C. 128 Wharton Road - Bristol, PA - 19007
uted By:	Customer Service: Toll Free: 877.877.9680 E-mail: sales@bioplastmfg.com
Distributed	TISO Margueret

©2015 Bioplast Manufacturing, L.L.C.