


Manufacturer of Quality Dishes Since 1968



Rim Feature	Code	Nominal Size (mm)	Style	ISO Mark Target
Slippable	2916	60 x 15	Mono	No



Dishes/ Sleeve	Dishes/	Box/ Pallet	Box Dimensions			Вох	Вох	
	Вох		L	W	Н	Cubic	Weight	
20 500	F00	0 72	12.25 in	12.50 in	12.50 in	1.11 ft³	9.22 lbs	
	500		31.12 cm	31.75 cm	31.75 cm	0.03 m ³	4.18 kgs	

KORD-Valmar

kord-valmark.com

Bioplast Manufacturing, L.L.C. 128 Wharton Road - Bristol, PA - 19007

Customer Service:

Toll Free: 877.877.9680 E-mail: sales@bioplastmfg.com

2916 - 60 x 15 mm Slippable Petri Dish

- Smaller diameter uses less media "Media Saver"
- Ideal for water treatment and pharmaceutical membrane tests
- Clear viewing surface
- Clear top and bottom viewing surfaces
- Ideal for use in automation
- Disposable (single use) polystyrene petri dish for general microbiological use
- Molded from medical-grade virgin polystyrene
- Optically clear
- Radiation sterilized
- Contents sterile if sleeve is unopened
- Ventilation ribs allow for free air circulation and reduce condensation during incubation
- Consistent flatness allows for even media distribution
- Flared lid skirt and squared corners for easy one-hand operation, especially with gloves
- Certificates of Sterilization Processing available upon request
- Traceable by lot and SKU number
- Certificates of Compliance available upon request

*media not included





Specially designed "slippable" lid is ideal for use in automated equipment



Smaller diameter uses less media -"Media Saver"



Ventilation ribs allow for free air circulation and reduce condensation during incubation







Manufacturer of Quality Dishes Since 1968



Manufactured to the highest quality for consistent flatness – radiation sterilized Ventilation Ribs allow for full air circulation and reduce condensation during incubation


Flared lid Skirt and squared corners for easy one-hand operation, especially with gloves

Nominal	K-V	Style	Rim	Unique Feature	Dishes/	Dishes/ Case	Cases/	
Size (mm)	e (mm) Code Style		Feature		Sleeve		Skid	
60 x 15	2901	Mono	Stackable	Media Saver	20	500	72	
100 x 15	2800	Mono	Fully Stackable	Stable Stacking Design	25	500	30	
100 x 15	2900	Mono	Stackable	The Industry Standard	25	500	30	
100 x 15	2903	Ві	Stackable	Two Compartments	25	500	30	
100 x 15	2912	Tri	Stackable	Three Compartments	25	500	30	
100 x 15	2913	Quad	Stackable	Four Compartments	25	500	30	
100 x 25	2906	Mono	Stackable	Extra Deep Dish	17	340	30	



Automation Petri Dishes (Slippable / No Rim Lid Designs)									
Nominal Size (mm)	K-V Code	Style	Rim Feature	Unique Feature	ISO Mark Target	Dishes/ Sleeve	Dishes/ Case	Cases/ Skid	
60 x 15	2916	Mono	Slippable	Media Saver	No	20	500	72	
95 x 15	2914	Mono	Slippable	Media Saver	No	25	500	30	
100 x 10	2904	Mono	Slippable	Space Saver	No	30	750	24	
100 x 15	1920	Mono	No Rim	Automated Lid	No	25	480	30	
100 x 15	2910	Mono	Slippable	The Industry Standard	Yes	25	500	30	
100 x 15	2917	Mono	Slippable	Clear Viewing Surfaces	No	25	500	30	
100 x 15	2920	Mono	No Rim	Automation Lid	No	20	500	30	
100 x 15	2911	Ві	Slippable	Two Compartments	Yes	25	500	30	
100 x 20	2905	Mono	Slippable	Deep Dish, 20mm	Yes	20	400	30	
100 x 25	2909	Mono	Slippable	Extra Deep for Seed Germination Tests	No	17	340	30	
85 x 13	2950	Mono	Fully Stackable/Slippable	New 85 x 13 Ultra Plate	Yes	30	600	30	
150 x 15	2902	Mono	No Rim	Concentric Bottom Ring Ensures Even Media Depth	No	10	100	40	





Contact Us!

Bioplast Manufacturing, L.L.C.

128 Wharton Road - Bristol, PA - 19007

Customer Service:

Toll Free: 877.877.9680

E-mail: sales@bioplastmfg.com

