

# Defender™ D500M

Mechanical Floor Beam Scale





### **Economical and Simple for General Weighing Needs**

The OHAUS Defender<sup>™</sup> D500M is a portable mechanical industrial scale designed to meet a variety of applications where simplicity and economy are a must. The D500M features an easy-to-read, dual-side flat beam with a shadow-proof poise, providing quick reading under most lighting conditions. A positive beam lock secures the indicator when the scale is not in use. Counterpoise weights allow for scale capacities of up to 1,000 pounds. The D500M scale utilizes a rugged fabricated steel design and is equipped with durable phenolic wheels for easy maneuverability.

#### Standard Features Include:


- Versatile Adaptable to in-plant, receiving, commercial and general weighing applications
- Durable Constructed of fabricated carbon steel and designed for accuracy and long life
- Maneuverable Lightweight and easy to maneuver with low resistance, wear-resistant phenolic wheels
- Legal for Trade The D500M is NTEP certified as a legal for trade device in the US

## **Defender**<sup>TM</sup> **D500M** Mechanical Beam Scale

#### **Specifications**

Model	D500M
Capacity x Readability	1,000 lb x 0.5 lb
NTEP Certified Resolution	1:2,000
Platform size	19.25 x 28.13 in
Beam	Die-cast dual-side flat beam
Base	Fabricated carbon steel
Poise	Shadow-proof positive beam lock
Wheels	4 x 5 in diameter, phenolic
Counterpoise Weights	Total 1,000 lb: 2 x 100 lb, 1 x 200 lb, 1 x 500 lb
Net Weight	113.5 lbs
Shipping Weight	124 lbs
Shipping Dimensions (W x D x H)	50.0 x 31.5 x 7.3 in

#### **Outline Dimensions**



#### **Approvals**

NTEP CC No. 08-070

Quality: OHAUS quality management system registered to ISO 9001:2008

#### **OHAUS CORPORATION**

\* 7 Campus Drive Suite 310 Parsippany, NJ 07054 USA

Tel: 800.672.7722 973.377.9000 Fax: 973.944.7177

#### www.ohaus.com

With offices throughout Europe, Asia, and Latin America

#### \* ISO 9001

Registered Quality Management System

