

5 Liter General Purpose Water Bath



WB05

Key Specifications

Description	Gen. Purpose Water Bath, Digital, 5 L
Display	3.25" LCD
Working Temperature Range °C	Ambient +5° to 100°
Temperature Stability °C	±0.1°
Temperature Calibration Capability	1-point
Working Access (L x W x D) (cm)	12.7 x 27.4 x 15.2
Flammability Class (DIN 12876-1)	I (NFL)

Key Features

Common Specifications

Working Temperature Range °F	Ambient +10° to 212°
Working Temperature Range °C	Ambient +5° to 100°
Reservoir Capacity (gallons)	1.32
Reservoir Capacity (liters)	5
Reservoir/Tank Material	Stainless Steel

Reservoir Cover	Hinged Gable
Working Access (L x W x D) (inches)	5 x 10.8 x 6
Working Access (L x W x D) (cm)	12.7 x 27.4 x 15.2
Temperature Stability °F	±0.2°
Temperature Stability °C	±0.1°
Display	3.25" LCD
Flammability Class (DIN 12876-1)	I (NFL)
Over-Temperature Protection / Failsafe Heater Control	Yes
Maximum Ambient Temperature °F	104°
Maximum Ambient Temperature °C	40°
Overall Dimensions (L x W x H) (inches)	12 x 14.5 x 13.5
Overall Dimensions (L x W x H) (cm)	30.5 x 36.8 x 34.5
Shipping Weight (pounds)	14.0
Shipping Weight (kilograms)	6.4
Catalog Page Number	86-87

50 Hz Only

Part Number	WB05A12E
Included Hardware	Sample Tray
Heater Wattage	360
Electrical Requirements (VAC/Hz/Ph/A)	240/50/1/1.6
Regulatory Approvals	CE

© 2015 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

Source URL (modified on 10/09/2017 - 13:13): <https://www.polyscience.com/general-purpose-water-baths/5-liter-general-purpose-water-bath>