

20 Liter General Purpose Water Bath



WB20

Key Specifications

Description	Gen. Purpose Water Bath, Digital, 20 L
Display	3.25" LCD
Working Temperature Range °C	Ambient +5° to 100°
Temperature Stability °C	±0.1°
Temperature Calibration Capability	1-point
Working Access (L x W x D) (cm)	24.1 x 43.2 x 15.2
Flammability Class (DIN 12876-1)	I (NFL)

Key Features

Common Specifications

Working Temperature Range °F	Ambient +10° to 212°
Working Temperature Range °C	Ambient +5° to 100°
Reservoir Capacity (gallons)	5.28
Reservoir Capacity (liters)	20
Reservoir/Tank Material	Stainless Steel

Reservoir Cover	Hinged Gable
Working Access (L x W x D) (inches)	9.5 x 17 x 6
Working Access (L x W x D) (cm)	24.1 x 43.2 x 15.2
Temperature Stability °F	±0.2°
Temperature Stability °C	±0.1°
Display	3.25" LCD
Flammability Class (DIN 12876-1)	I (NFL)
Over-Temperature Protection / Failsafe Heater Control	Yes
Reservoir Drain	Yes
Maximum Ambient Temperature °F	104°
Maximum Ambient Temperature °C	40°
Overall Dimensions (L x W x H) (inches)	17.5 x 22.5 x 15.5
Overall Dimensions (L x W x H) (cm)	44.5 x 57.2 x 39.3
Shipping Weight (pounds)	25.0
Shipping Weight (kilograms)	11.3
Catalog Page Number	86-87

50 Hz Only

Part Number	WB20A12E
Included Hardware	Sample Tray
Heater Wattage	1400
Electrical Requirements (VAC/Hz/Ph/A)	240/50/1/6
Regulatory Approvals	CE

© 2015 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

Source URL (modified on 10/09/2017 - 13:14): <https://www.polyscience.com/general-purpose-water-baths/20-liter-general-purpose-water-bath>