

Lab Ovens

LT - Low-Temp Series Ovens from Quincy Lab, Inc.





FEATURING

■ Digital Microprocessor Control
 ■ Air Forced Gravity Convection
 ■ Durable Double-wall
 Construction
 ■ Scratch-resistant Hammer Finish
 ■ Corrosion-resistant Aluminized Interior
 ■ Full 18-month Factory Warranty
 ■ Proven Reliability

Advanced Design

Cabinets are as attractive as they are durable. Exteriors are painted light gray and have a hard, scratch-resistant hammer finish. Doors open with high-impact thermoplastic handles. Cabinets have heavy steel double-wall construction. Work space is insulated from the outer cabinet with one inch of high-density mineral wool, and interiors are made of corrosion-resistant aluminized steel.

Controls

The digital control combines the features of the analog model but offers the ease of temperature setting and the stability of a full PID microprocessor that accurately maintains settings within +/- 1.0° C, even in varying ambient or power supply

conditions. A temperature-tracking feature stores temperature deviation from set point. This feature helps to confirm stability or indicate any control malfunction or power loss throughout a process period. The digital controller features large LED's that continuously display process temperature as well as a setting lock mode that provides protection against accidental or inadvertent adjustment.

Heating Elements

Energy-efficient, low-watt density incoloy sheathed elements are engineered into a compact design for quick run-up and recovery times. Temperature uniformity is greatly improved by a perforated heat shield which absorbs radiant heat and distributes it more evenly.


GCE-LT and AFE-LT Models

Quincy Lab uses only quality UL and CSA recognized components in all ovens. The 'GCE-LT' series gravity convection ovens are designed to meet the laboratory needs of industry, research organizations, and schools. Well-crafted and versatile, they are used for part drying, baking, curing, sterilizing, evaporating, heat treating, annealing, and testing. The 'AFE-LT' models offer forced-air circulation allowing for uniform distribution of heat throughout the chamber, which also facilitates dehydration and evaporating. The GCE-LT and AFE-LT series ovens have temperature ranges of 210°F/99°C, and 225°F/107°C respectively, and are competitively priced and offer exceptional value and reliability.


Quincy Lab, In. has been a mainline manufacturer of laboratory ovens and incubators for more than 40 years. We are dedicated to product value, customer satisfaction, and ongoing product support.

GENERAL SPECIFICATIONS	MODEL 10GCE-LT	MODEL 20GCE-LT	MODEL 30GCE-LT	MODEL 40GCE-LT	MODEL 10AFE-LT	MODEL 20AFE-LT	MODEL 30AFE-LT	MODEL 40AFE-LT
INTERIOR DIMENSIONS INCHES W x H x D (CM) W x H x D	12x10x10 30.5x25.4x25.4	13x13x13 33x33x33	18x15.6x12 45.7x39.6x30.5	18x21.8x14 45.7x55.4x35.6	12x8.25x10 30.5x21x25.4	13x11x13 33x28x33	18x14.2x12 45.7x36x30.5	18x19.8x14 45.7x50.3x35.6
EXTERIOR DIMENSIONS INCHES W x H x D (CM) W x H x D	14x17.5X12.3 35.6x44.5x31.2	15x21.5x15.3 38x54.6x38.9	20x25.5x14.3 50.8x64.8x36.2	20x31.5X16.3 50.8x80x41.4	14x20.5X12.3 35.6x52x31.2	15x25x15.3 38x63.5x38.9	20x29x14.3 50.8x73.7x36.2	20x35X16.3 50.8x89x41.4
CAPACITY CUBIC FEET (LITERS)	0.7 19.8	1.27 36	2.0 56.6	3.0 85	0.6 17	1.14 32.3	1.83 51.8	2.86 81
TEMPERATURE RANGE AMBIENT +15F TO F / C	210/99	210/99	210/99	210/99	225/107	225/107	225/107	225/107
SHELVES (1" CENTERS) MAXIMUM PER UNIT MAXIMUM POUNDS / SHELF	10 35	13 35	16 35	22 35	8 35	11 35	14 35	20 35
ELECTRICAL* VOLTS/AMPS WATTS PLUG/NEMA	115 / 2.6 360 5-15P	115 / 2.6 360 5-15P	115 / 5.2 720 5-15P	115 / 5.2 720 5-15P	115 / 3.65 440 5-15P	115 / 3.65 440 5-15P	115 / 4.8 800 5-15P	115 / 6.8 800 5-15P
WEIGHT SHIPPING STAND ALONE	42.6 31.4	57 43.5	76.2 57.2	91.5 72.5	50.4 37.2	65.7 58.7	84.2 64.2	98.5 80.8


^{*} Standard models voltage only, 230 voltage available.



'GCE-LT' MODEL AIR FLOW



'AFE-LT' MODEL AIR FLOW



**PERFORMANCE CHARACTERISTICS	MODEL 10GCE-LT	MODEL 20GCE-LT		MODEL 40GCE-LT	MODEL 10AFE-LT	MODEL 20AFE-LT	MODEL 30AFE-LT	MODEL 40AFE-LT
CONTROL STABILITY @ 75°C	+/- 0.8°C	+/- 0.8°C	+/- 1.0°C	+/- 1.0°C	+/- 0.8°C	+/- 0.8°C	+/- 1.0°C	+/- 1.0°C
TIME TO TEMPERATURE AMBIENT TO MAX	27 Min.	30 Min.	28 Min.	34 Min	36 Min.	40 Min.	27 Min.	40 Min.
RECOVERY @200° F DOOR OPEN 15 SEC. DOOR OPEN 30 SEC.	4 Min. 8 Min.	5.5 Min. 8 Min.	5 Min. 8 Min.	8 Min. 14 Min.	3 Min. 4 Min.	8 Min. 12 Min.	4 Min. 6 Min.	7 Min. 9 Min.

^{**} PERFORMANCE CHARACTERISTICS FOR STANDARD VOLTAGE MODELS, ALTERNATE VOLTAGE MODELS MAY VARY. ALL TESTS CONDUCTED UNDER CONTROLLED LABORATORY CONDITIONS.

