

Product Specification

Alkaline-Iodide-Azide for Dissolved Oxygen (DO) Analysis

Lot Number: SAMPLE

Product Number: 541

Manufacture Date: N/A

Expiration Date: N/A

This solution is intended for use with supersaturated samples. White turbidity due to Sodium Carbonate may form during testing, but this will not affect the results of the test.

Name	CAS#	Grade
Sodium Iodide	7681-82-5	ACS
Water	7732-18-5	ACS/ASTM/USP/EP
Sodium Hydroxide	1310-73-2	ACS
Sodium Azide	26628-22-8	Reagent

Test	Specification	Result
Appearance	Colorless liquid	N/A
Free Iodine	To Pass Test	N/A

Specification	Reference
Alkali-Iodide-Azide Reagent, for supersaturated samples	APHA (4500-O D)
Alkali-Iodide-Azide Reagent, for supersaturated samples	APHA (4500-O E)
Alkali-Iodide-Azide Reagent, for supersaturated samples	APHA (4500-O F)
Alkali-Iodide-Azide Reagent, for supersaturated samples	APHA (4500-O C)

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
541-32	1 L natural poly	24 months
541-16	500 mL natural poly	24 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."