

Product Specification

Reagent Alcohol, 95% (v/v)

Lot Number: SAMPLE**Product Number:** R2546950**Manufacture Date:** N/A**Expiration Date:** N/A

Name	CAS#	Grade
Ethyl Alcohol	64-17-5	Reagent
Water	7732-18-5	Reagent
Methyl Alcohol	67-56-1	Reagent
Isopropyl Alcohol	67-63-0	Reagent

Test	Specification	Result
Appearance	Colorless liquid	N/A
Specific Gravity	0.810-0.812 at 20°C	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
R2546950-4A	4 L natural poly	36 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."