

Product Specification

Ferric Chloride TS, 9% (w/v)

Lot Number: SAMPLE**Product Number:** 3109**Manufacture Date:** N/A**Expiration Date:** N/A

This solution may become cloudy or hazy due to the absence of acid. Keep the container tightly closed.

Name	CAS#	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Ferric Chloride Hexahydrate	10025-77-1	ACS

Test	Specification	Result	NIST SRM#
Appearance	Bright orange/brown liquid	N/A	
Assay (vs. Sodium Thiosulfate/Starch)	8.9-9.1 % (w/v)	N/A	136

Specification	Reference
Ferric Chloride TS	USP (N/A)
Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with ϵ thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.	

Part Number	Size / Package Type	Shelf Life (Unopened Container)
3109-4	120 mL amber glass	6 months
3109-16	500 mL amber glass	6 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."