

Product Specification

Hydrogen Peroxide, 10% (w/w), Reagent Grade Stabilized with Acetanilide

Lot Number: SAMPLE

Product Number: R3821310

Manufacture Date: N/A

Expiration Date: N/A

Name	CAS#	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Hydrogen Peroxide	7722-84-1	ACS
Acetanilide	103-84-4	Reagent

Test	Specification	Result	NIST SRM#
Appearance	Colorless liquid	N/A	
Assay (vs. Potassium Permanganate)	9.9-10.1 % (w/w)	N/A	40

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
R3821310-1BV	1 L vented opaque poly	18 months
R3821310-20F	20 L Cubitainer®	18 months
R3821310-4BV	4 L vented black poly	18 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)