

Product Specification

Stannous Chloride Reagent I for Phosphorus (Phosphate) Analysis (without extraction)

Lot Number: SAMPLE

Product Number: 7997

Manufacture Date: N/A

Expiration Date: N/A

Name	CAS#	Grade
Glycerin	56-81-5	ACS
Stannous Chloride Dihydrate	10025-69-1	ACS

Test	Specification	Result
Appearance	Colorless to light tan liquid	N/A

Specification	Reference
Stannous Chloride Reagent I	APHA (4500-P D)

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
7997-4	120 mL natural poly	18 months
7997-32	1 L natural poly	18 months
7997-16	500 mL natural poly	18 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)