

Product Specification

**Triton® X-100
Wetting Agent**
Lot Number: SAMPLE

Product Number: 8698.5

Date Packaged: N/A

Expiration Date: N/A

Triton® is a registered trademark of Union Carbide Corp.
Typical Values (not specifications):
Color (APHA): max. 100
Specific gravity: 1.065 at 25°C
Viscosity (Brookfield): 240 cps at 25°C
Cloud point (1% solution): 63-69°C

Contains
 $(\text{CH}_3)_3\text{CCH}_2\text{C}(\text{CH}_3)_2\text{C}_6\text{H}_4(\text{OCH}_2\text{CH}_2)_{9,10}\text{OH}$ CAS No. 9002-93-1 MW 624.83

Test	Specification	Result
Appearance	Colorless liquid	N/A

Specification	Reference
Standard Isotonic Sodium Chloride Solution for Immuno-hematologic Testing	NCCLS (ASI/BB-1)

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
8698.5-32	1 L natural poly	36 months
8698.5-5	20 L Cubitainer®	36 months
8698.5-1	4 L natural poly	36 months
8698.5-16	500 mL natural poly	36 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."