

Product Specification

Sulfuric Acid, Technical, 66° Bé

Lot Number: SAMPLE**Product Number:** RABS0050**Date Packaged:** N/A**Expiration Date:** N/A**Contains** H_2SO_4 CAS No. 7664-93-9 MW 98.08

Test	Specification	Result
Assay	min 93.0 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RABS0050-2.5D	2.5 L safety-coated glass	36 months
RABS0050-4x2.5D	4 x 2.5 L safety-coated glass	36 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."