

Product Specification

Ammonium Thiocyanate, ACS Reagent Grade

Lot Number: SAMPLE**Product Number:** RDCA0550**Date Packaged:** N/A**Expiration Date:** N/AChemical formula: NH_4SCN
Formula weight: 76.12**Contains** NH_4SCN CAS No. 1762-95-4 MW 76.12

Test	Specification	Result
Assay	min 97.5 %	N/A
Chloride (Cl)	max 0.005 %	N/A
Heavy Metals (as Pb)	max 5 ppm	N/A
Insoluble Matter	max 0.005 %	N/A
Iodine-consuming substances	max 0.004 meq/g	N/A
Iron (Fe)	max 3 ppm	N/A
pH of 5% Solution	4.5-6.0	N/A
Residue after Ignition	max 0.025 %	N/A
Sulfate (SO_4)	max 0.005 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCA0550-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."