

Product Specification

Dichloromethane, ACS Reagent Grade (Methylene Chloride)

Lot Number: SAMPLE

Product Number: RDCD0100

Date Packaged: N/A

Expiration Date: N/A

Contains
CH ₂ Cl ₂ CAS No. 75-09-2 MW 84.93

Test	Specification	Result
Appearance	Colorless liquid	N/A
Assay	min 99.5 %	N/A
Color (APHA)	max 10	N/A
Free Halogens	To Pass Test	N/A
Residue after evaporation (Non-Volatile Matter)	max 0.002 %	N/A
Titration Acid	max 0.0003 meq/g	N/A
Water (H ₂ O)	max 0.02 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCD0100-1C	1 L amber glass	36 months
RDCD0100-4C	4 L amber glass	36 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."