

Product Specification

p-Dimethylaminobenzaldehyde, ACS Reagent Grade

Lot Number: SAMPLE**Product Number:** RDCD0500**Date Packaged:** N/A**Expiration Date:** N/A**Contains** $C_9H_{11}NO$ CAS No. 100-10-7 MW 149.19

Test	Specification	Result
Color (APHA) of alcohol solution	max 60	N/A
Color of Hydrochloric Acid solution	To Pass Test	N/A
Melting Point	73-75 °C	N/A
Residue after Ignition	max 0.1 %	N/A
Solubility in Alcohol	To Pass Test	N/A
Solubility in Hydrochloric Acid	To Pass Test	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCD0500-100B1	100 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."