

Product Specification

Lead Acetate Trihydrate, ACS Reagent Grade

Lot Number: SAMPLE

Product Number: RDCL0150

Date Packaged: N/A

Expiration Date: N/A

Contains $(\text{CH}_3\text{COO})_2\text{Pb} \cdot 3\text{H}_2\text{O}$ CAS No. 6080-56-4 MW 379.33

Test	Specification	Result
Assay	99.0-103.0 %	N/A
Calcium (Ca)	max 0.005 %	N/A
Chloride (Cl)	max 5 ppm	N/A
Copper (Cu)	max 0.002 %	N/A
Insoluble Matter	max 0.01 %	N/A
Iron (Fe)	max 0.001 %	N/A
Nitrate and Nitrite (as NO_3)	max 0.005 %	N/A
Potassium (K)	max 0.005 %	N/A
Sodium (Na)	max 0.01 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCL0150-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)