

Product Specification

1,10-Phenanthroline Monohydrate, ACS Reagent Grade

Lot Number: SAMPLE**Product Number:** RDCP0200**Date Packaged:** N/A**Expiration Date:** N/A**Contains** $C_{12}H_8N_2 \cdot H_2O$ CAS No. 5144-89-8 MW 198.22

Test	Specification	Result
Suitability as Redox Indicator	To Pass Test	N/A
Suitability for determining Iron	To Pass Test	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCP0200-5F1	5 g pouch	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."