

Product Specification

Picric Acid, Reagent Grade

For safety in transportation, Picric Acid is usually mixed with a minimum of 30% water.

Lot Number: SAMPLE

Product Number: RDCP0300

Date Packaged: N/A

Expiration Date: N/A

Contains

$(\text{NO}_2)_3\text{C}_6\text{H}_2\text{OH}$ CAS No. 88-89-1 MW 229.10

Test	Specification	Result
Identification	To Pass Test	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCP0300-10C1	10 g amber glass	84 months
RDCP0300-125A6	125 g poly natural	84 months
RDCP0300-25C1	25 g amber glass	84 months
RDCP0300-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."