

Product Specification

Potassium Antimony Tartrate Trihydrate, ACS Reagent Grade

Lot Number: SAMPLE

Product Number: RDCP0320

Date Packaged: N/A

Expiration Date: N/A

Contains $K_2(C_4H_4O_6Sb)_2 \cdot 3H_2O$ CAS No. 11071-15-1 MW 671.90

Test	Specification	Result
Arsenic (As)	max 0.015 %	N/A
Assay	99.0-103.0 %	N/A
Loss on drying	max 2.7 %	N/A
Titration acid or base	max 0.020 meq/g	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCP0320-100B1	100 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."