

Product Specification

Potassium Hydroxide, ACS Reagent Grade

Lot Number: SAMPLE**Product Number:** RDCP0500**Date Packaged:** N/A**Expiration Date:** N/AChemical formula: KOH
Formula weight: 56.11**Contains**

KOH CAS No. 1310-58-3 MW 56.11

Test	Specification	Result
Assay	min 85 %	N/A
Calcium (Ca)	max 0.005 %	N/A
Chloride (Cl)	max 0.01 %	N/A
Heavy Metals (as Ag)	max 0.001 %	N/A
Iron (Fe)	max 0.001 %	N/A
Magnesium (Mg)	max 0.002 %	N/A
Nickel (Ni)	max 0.001 %	N/A
Nitrogen Compounds (as N)	max 0.001 %	N/A
Phosphate (PO ₄)	max 5 ppm	N/A
Potassium Carbonate (K ₂ CO ₃)	max 2.0 %	N/A
Sodium (Na)	max 0.05 %	N/A
Sulfate (SO ₄)	max 0.003 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCP0500-1000B1	1 kg white poly	84 months
RDCP0500-100B1	100 g poly	84 months
RDCP0500-3A5	3 kg poly natural	84 months
RDCP0500-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."