

Product Specification

Sodium Sulfite, Anhydrous, ACS Reagent Grade

Lot Number: SAMPLE

Product Number: RDSCS0870

Date Packaged: N/A

Expiration Date: N/A

Contains Na_2SO_3 CAS No. 7757-83-7 MW 126.04

Test	Specification	Result
Assay	min 98.0 %	N/A
Chloride (Cl)	max 0.02 %	N/A
Free acid	To Pass Test	N/A
Heavy Metals (as Pb)	max 0.001 %	N/A
Insoluble Matter	max 0.005 %	N/A
Iron (Fe)	max 0.001 %	N/A
Titration Free Base	max 0.03 meq/g	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDSCS0870-12F4	12 kg fiber drum	84 months
RDSCS0870-2.5A5	2.5 kg poly natural	84 months
RDSCS0870-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."