

Product Specification

Sodium Thiocyanate, ACS Reagent Grade

Lot Number: SAMPLE

Product Number: RDCS0910

Date Packaged: N/A

Expiration Date: N/A

Contains

NaSCN CAS No. 540-72-7 MW 81.07

Test	Specification	Result
Ammonium (NH ₄)	max 0.002 %	N/A
Assay	min 98.0 %	N/A
Carbonate (as Na ₂ CO ₃)	max 0.2 %	N/A
Chloride (Cl)	max 0.01 %	N/A
Heavy Metals (as Pb)	max 5 ppm	N/A
Insoluble Matter	max 0.005 %	N/A
Iron (Fe)	max 2 ppm	N/A
Sulfate (SO ₄)	max 0.01 %	N/A
Sulfide (S)	max 0.001 %	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RDCS0910-100B1	100 g poly	84 months
RDCS0910-500B1	500 g poly	84 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)