

Product Specification

Kerosene - Odorless**Lot Number:** SAMPLE**Product Number:** RSOK0100**Date Packaged:** N/A**Expiration Date:** N/A**Contains** $C_{6-16}H_{12-34}$ CAS No. 8008-20-6

Test	Specification	Result
Appearance	Colorless liquid	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RSOK0100-1C	1 L amber glass	60 months
RSOK0100-4C	4 L amber glass	60 months
RSOK0100-500C	500 mL amber glass	60 months
RSOK0100-55E	55 gal drum	60 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."