

Product Specification

Morpholine, ACS Reagent Grade

Lot Number: SAMPLE

Product Number: RSOM0200

Date Packaged: N/A

Expiration Date: N/A

Contains C_4H_9NO CAS No. 110-91-8 MW 87.12

Test	Specification	Result
Assay	min 99.0 %	N/A
Boiling Range	126 - 130 °C	N/A
Color (APHA)	max 15	N/A

Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document traceability and production history for each lot manufactured.

Part Number	Size / Package Type	Shelf Life (Unopened Container)
RSOM0200-250C	250 mL amber glass	60 months
RSOM0200-500C	500 mL amber glass	60 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)